

Plinth 2000

Model 503 3 Section Couch

Overview:

Very popular and versatile three-section couch, firm and stable providing the ideal base for general treatment techniques through to heavy manipulation. Featuring smooth anti-entrapment gas operated back & foot rest, lifetime frame and 5 year parts guarantee, this couch forms the backbone of many clinics. A wide range of options are available to expand the use of the couch.




Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Hoist clearance 130mm
- Anti-entrapment back & foot rest

Overall Dimensions

Length: 183cm (72") approx
Width: 63cm (25") approx
Height: 47cm - 97cm approx
(19" - 38") approx

Upholstery Dimensions



Recommended Options:



Double Foam




Battery Back-up



Folding Cotside



Motor Back



Breathing Hole



Paper Roll Holder


NHS Approved Supplier



UK Made

t: +44 (0)1449 767887
e: sales@plinth2000.com
www.plinth2000.com

Operational Parameters



Model 503E - Electric version

Model 503H - Hydraulic version

For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



Plinth 2000 Ltd.
Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk, IP14 5PP

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com
www.plinth2000.com

Plinth | 2000