

Mobile Dialysis
Chair
Platinum Range



PLATINUM RANGE

MODEL 94DY

Developed specifically for busy renal units, this dual column chair features vertical lifting to minimise space requirements and features simple push-button electric operation of all movement parameters. These features include electric profiling of height, backrest and foot section angles, foot support board position, Trendelenburg tilt, auto CPR recovery position and return to chair position. These benefits provide ideal solutions for handling high risk patients for procedures such as dialysis, periton dialysis, plasmapheresis diagnosis, oncology and chemotherapy.

CPR

In the event of an emergency, the dialysis chair has an electric CPR function to reposition the chair as quickly as possible for the procedure.



Seat Tilt

By combining the electric operation of height, tilt and backrest angles, various seated and semi-seated or reclined positions can be achieved for the benefit of the patient or for the required position for various treatments and procedures. The fully upholstered armrests provide patient support whilst using the chair. A matching removable head support cushion can also be added to provide further comfort.



Trendelenburg Tilt

The Trendelenburg tilt facility is electrically operated with a range up to -20 degrees.



Transport Mode

The couch has a large central locking wheel system with steering facility for transport procedures and can be used in both the seated or trolley modes.

Adjustable Armsrests

The fully upholstered armrests provide patient support in the seated or semi-seated position and automatically adjust to the arm support position as the backrest angle is electronically positioned. Each armrest is width and angle adjustable to suit individual user requirements. A matching removable head support cushion can also be added to provide additional patient comfort.



MODEL 94DY

Product Features

Dual Column Lift	Vertical lifting to minimise space
Electric Operation	Controls height, tilt, backrest & footrest angles, CPR and Auto-Return
Safe Working Load	250kgs (550lbs or 39 stone) approx.
Height Range	52cm to 102cm (20" to 40") approx.
Backrest Angle	Horizontal to +80 degrees
Foot Section Angle	Electric operation from +30 degrees to -35 degrees
Trendelenburg Tilt	Electric operation from 0 to -20 degrees
Wheel System	Large 125mm twin castors with central locking and steering facility
Footboard	Electrically operated - extends 28cm (11")
Arm Supports	Width adjustable with auto-angular adjustment
Battery Back-Up	Providing emergency power during mains power failure
CPR Facility	Electrically operated via hand switch control module
Hoist Clearance	Clearance beneath frame for mobile hoist - 13cm (5")
Upholstery	Deep sculptured pressure relieving design with antimicrobial vinyl
Dimensions (LxWxH)	198cm x 60cm x 52-102cm (78" x 24" x 20-40") approx.
Warranty	Lifetime frame and 5 year parts guarantee

Platinum Range Products:

- Dialysis
- Podiatry
- Oncology
- Gynaecology
- Cardio
- Bariatric
- Lymphoedema
- Ophthalmology
- Leg Ulcer
- Ultrasound

Plinth | Medical

Plinth Medical Ltd.
Barric Lane,
Occold, Eye, Suffolk,
IP23 7PX

+44 (0)1449 767887
sales@plinthmedical.com

www.plinthmedical.com

NHS Approved Supplier

