

**Plinth** | Medical

Mobile  
Cardio Station  
Platinum Range



# PLATINUM RANGE MODEL 94MCS

Developed for busy Outpatient Departments such as Cardiology, this dual column couch features vertical lifting to minimise space requirements and provides all the benefits of a mobile transport chair or patient trolley, with all-electric profiling in both modes. The couch has a safe working load of 250kgs and additional features such as electric CPR mode, Trendelenburg tilt, battery back-up and auto-return to chair set-up position.

## Trendelenburg Tilt

The Trendelenburg tilt facility is electrically operated with a range up to -20 degrees. The fold-down padded side support rails provide additional security for this procedure.



## Seat Tilt

The combination of electric height, tilt and backrest operation allows various seated and semi-seated or reclined positions to be achieved. Such positions can be for the benefit of the patient or for achieving the required position for the treatment or procedure that they are undergoing.



## Transport Handles

The transport handles are positioned beneath the backrest section and are turned inwards for storage when not in use. They are quickly and easily extended and locked into position when required for transport purposes. The backrest also provides a storage pocket for patient notes or belongings.

## Mobile Chair Mode

The couch is electrically positioned for seated procedures and the transport handles can be positioned and set into the required position for patient transport. The shortened side rails, with their arm rests and side pads add both security and comfort for the individual user. In addition, the backrest and foot sections can be adjusted to various angles to improve patient comfort.



## Trolley Mode

The Cardio station can be quickly configured into the transport trolley mode for patient positioning during rest or recovery periods or for transport from clinic to clinic or to other locations around the hospital.



# MODEL 94MCS

## Product Features

Dual Column Lift	Vertical lifting to minimise space
Electric Operation	Controls height, tilt, backrest/footrest angles, CPR and Auto-Return
Safe Working Load	250kgs (550lbs or 39 stone)
Height Range	52cm to 102cm (20" to 40") approx.
Backrest Angle	Horizontal to +80 degrees
Foot Section Angle	Electric operation from +10 degrees to -90 degrees
Trendelenburg Tilt	Electric operation from 0 to -20 degrees
Wheel System	Large 125mm twin castors with central locking and steering facility
Side Rails	Extending side rails (seat to trolley) with side pads and armrests
Push Handles	Fully extending and lockable; stowaway storage facility
Battery Back-Up	Providing emergency power during mains power failure
CPR Facility	Electrically operated via hand switch control module
Hoist Clearance	Clearance beneath frame for mobile hoist - 13cm (5")
Upholstery	Pressure relieving foam design with antimicrobial vinyl - 75mm deep
Dimensions (LxWxH)	193cm x 79cm x 52-102cm (LxWxH) (77" x 31" x 20-40") approx.
Warranty	Lifetime frame and 5 year parts guarantee

### Platinum Range Products:

- Dialysis
- Podiatry
- Oncology
- Gynaecology
- Cardio
- Bariatric
- Lymphoedema
- Ophthalmology
- Leg Ulcer
- Ultrasound

**Plinth** | Medical

Plinth Medical Ltd.  
Barric Lane,  
Occold, Eye, Suffolk,  
IP23 7PX

+44 (0)1449 767887  
sales@plinthmedical.com

[www.plinthmedical.com](http://www.plinthmedical.com)

NHS Approved Supplier

