



### Overview:

This versatile and comfortable examination and treatment chair has been specifically designed for use in Ophthalmology and ENT Departments. The chair features electric operation of height, backrest and foot section angles. The backrest also provides a moulded head support with height adjustment for patient positioning. The backrest can be lowered to  $-15^{\circ}$  for patient recovery and the foot section lowers to  $-90^{\circ}$ . The chair has fold-down padded side supports for patient security and the large braked castors enable the chair to be easily moved and positioned within the department or clinic.

### Benefits


- Lifetime frame & 5 year parts guarantee - see website for further details
- Electric operation of height, backrest and foot section
- The backrest has a moulded patient head support with height adjustment
- The backrest is adjustable up to  $85^{\circ}$  for sitting procedures and lowers to  $-15^{\circ}$  for patient recovery
- Foot section adjustment is 0 to  $-90^{\circ}$  for patient comfort and access
- Fold-down padded side supports for patient comfort and security
- Safe working load of 250kgs (550lb / 39 stone)
- 125mm braked castors for easy transport and positioning of the chair
- Low patient transfer height of 52cm (20")
- Wide choice of options, width and accessories to meet user requirements




### Overall Dimensions

Length: 158cm (62") approx  
 Width: 63cm (25") approx  
 Height: 52cm - 101cm approx  
 (20" - 40") approx

### Upholstery Dimensions




### Options & Accessories:



Double Foam



Retractable Castors



Battery Back-up



Classic Operators Chair




Foot Switch Operation

NHS Approved Supplier



t: +44 (0)1449 767887  
 e: sales@plinthmedical.com  
 www.plinthmedical.com

## Operational Parameters



Model 930 - Electric Version

## Standard Features

- Designed and manufactured in the UK
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- 125mm toe braked castors
- Hoist clearance 130mm
- Motorised back and leg rest
- Adjustable locking head rest
- Fold down padded support rails

## Technical Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

## Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

## Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 - CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond