

Overview:

An extremely strong and rigid 2 section couch designed for various examination and manipulation procedures and suitable for constant heavy use. This model features a full length head section with breathing hole and plug which is ideal for supporting a patient in a sitting position. The head section provides both positive and negative positioning.



Model shown with the optional dual foot switch

Benefits

- Lifetime frame & 5 year parts guarantee - see website for further details
- Proven lifespan in excess of 16 years
- Suitable for heavy manipulations & high usage
- Low running & maintenance costs
- Customisable to suit user's needs
- Safe working load of 250kgs (550lb / 39 stone)
- Optional foot bar switch system provides all-round, hands free height operation
- Extendable foot section option increase overall couch length for taller patients
- Breathing hole and plug for prone procedures
- Head section can be lowered to -25° for head and neck procedures or for patient recovery

Overall Dimensions

Length: 195cm (77") approx
with foot extension 195 - 220cm (77" - 87") approx
Width: 63cm or 70cm
(25" or 28") approx
Height: 47cm - 97cm approx
(18" - 38") approx

Upholstery Dimensions



Options & Accessories:



Paper Roll Holder



Linked Castor Base



Foot Switch Operation



Foot Switch Extension



Dual Foot Switch Operation




Foot Bar Switch

NHS Approved Supplier



t: +44 (0)1449 767887
e: sales@plinthmedical.com
www.plinthmedical.com

Operational Parameters



Model 512E - Electric Version

Model 512H - Hydraulic Version

Standard Features

- Designed & manufactured in the UK
- Lifetime frame & 5 year parts guarantee
- Very strong & rigid frame
- High-performance bearings on all pivot points
- Safe working load & lift capacity 250kg
- Hydraulic or electric height operation
- Retractable wheel system
- Adjustable foot for uneven floors
- Anti-entrapment gas struts
- Head section adjustable from -25° to 80°
- Option of foot extension for taller patients
- Breathing hole and plug
- Choice of width 63cm or 70cm
- Hoist clearance 13cm

Technical Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 - CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Plinth Medical Ltd
Barric Lane
Occold, Eye
Suffolk, United Kingdom
IP23 7PX

t: +44 (0)1449767887
e: sales@plinthmedical.com
www.plinthmedical.com

Plinth | Medical