

Overview:

Designed for rehabilitation techniques, this model features a head section with breathing hole and plug, and height adjustable armrests, providing additional comfort for the patient during prone head and neck procedures. The leg section can be extended for the taller patient and has angular adjustment between -25° and +85°. Electric models are supplied as standard with hand switch control, with options of single or dual footswitch operation, or the hands-free foot bar switch system.



Model 515E with optional chrome foot bar switch system, linked castor base and foot extension.

Colour: Pillarbox Red


Benefits

- Lifetime frame & 5 year parts guarantee - see website for further details
- Proven lifespan in excess of 16 years
- Suitable for heavy manipulations & high usage
- Low running & maintenance costs
- Customisable to suit user's needs
- Optional foot bar switch system provides all-round, hands-free height operation
- Safe working load of 250kgs (550lb / 39 stone)
- Extendable foot section option increases overall length for taller patients
- Head section design and mid-section cut out enhances patient comfort for head and neck manipulation and mobilisation techniques

Overall Dimensions

Length: 186cm (72") approx
with foot extension 186 - 211cm (72" - 83") approx
Width: 63cm or 70cm (25" or 28") approx
Height: 47cm - 97cm approx
(18" - 38") approx

Upholstery Dimensions



Options & Accessories:



Paper Roll Holder



Foot Switch Operation



Dual Foot Switch Operation



Foot Bar Switch



Linked Castor Base




Foot Extension

NHS Approved Supplier



t: +44 (0)1449 767887
e: sales@plinthmedical.com
www.plinthmedical.com

Operational Parameters



Model 515E - Electric Version

Model 515H - Hydraulic Version

Standard Features

- Manufactured & designed in the UK
- Lifetime frame & 5 year parts guarantee
- Very strong & rigid frame
- High-performance bearings on all pivot points
- Safe working load & lift capacity 250kg
- Hydraulic or electric height operation
- Retractable wheel system
- Adjustable foot for uneven floors
- Anti-entrapment gas struts
- Extendable foot section
- Height-adjustable armrests
- Foot section adjustable from -25° to +85°
- Head section adjustable from -30° to +75°
- Cutaway mid-section for armrest access
- Choice of width 63cm or 70cm
- Hoist clearance 13cm

Technical Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 - CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

Plinth Medical Ltd
Barric Lane
Occold, Eye
Suffolk, United Kingdom
IP23 7PX

t: +44 (0)1449767887
e: sales@plinthmedical.com
www.plinthmedical.com

Plinth | Medical