



Overview:

A very strong and rigid couch designed for various manipulation and mobilisation procedures and mobilisation techniques. Suitable for constant and heavy usage, this couch will last for many years providing a trouble-free long life with minimal maintenance required. The short head section with breathing hole and plug is ideal for head and neck procedures.



Model shown with the optional foot extension


Benefits

- Lifetime frame & 5 year parts guarantee - see website for further details
- Proven lifespan in excess of 16 years
- Suitable for heavy manipulations & high usage
- Low running & maintenance costs
- Customisable to suit user's needs
- Optional foot bar switch system provides all-round, hands free height operation
- Safe working load of 250kgs (550lb / 39 stone)
- Extendable foot section option increase overall couch length for taller patients
- Head section can be lowered to -30° for head and neck procedures or for patient recovery

Overall Dimensions

Length: 186cm approx (73") approx
 with foot extension 186 - 211cm (73" - 83") approx
 Width: 63cm or 70cm (25" or 28") approx
 Height: 47cm - 97cm (18" - 38") approx

Upholstery Dimensions (cm)



Options & Accessories:



Paper Roll Holder



Linked Castor Base



Foot Switch Operation



Dual Foot Switch Operation



Foot Bar Switch



Foot Extension

Operational Parameters



Model 513SE - Electric Version

Model 513SH - Hydraulic Version

Standard Features

- Designed & manufactured in the UK
- Lifetime frame & 5 year parts guarantee
- Very strong & rigid frame
- High-performance bearings on all pivot points
- Safe working load & lift capacity 250kg
- Hydraulic or electric height operation
- Retractable wheel system
- Adjustable foot for uneven floors
- Anti-entrapment gas struts
- Short head section for head and neck procedures
- Head section adjustable from -30° to +75°
- Foot section adjustable from -25° to +85°
- Breathing hole and plug
- Choice of width: 63cm or 70cm
- Hoist clearance 13cm

Technical Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 - CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

Plinth Medical Ltd
Barric Lane
Occold, Eye
Suffolk, United Kingdom
IP23 7PX

t: +44 (0)1449767887

e: sales@plinthmedical.com

www.plinthmedical.com

Plinth | Medical