



### Overview:

The 50+ model provides an extra wide treatment surface of 102cm with a wider base frame for added stability and long life. This couch has electric operation of both height and backrest angle to assist clinicians working with bariatric and very heavy patients. The upholstery has the option of being upgraded with a pressure relieving membrane and multi-stretch vinyl to provide additional comfort for the patient.




### Benefits

- Lifetime frame & 5 year parts guarantee - see website for further details
- Extra large treatment surface of 189 x 102cm to support both patient and therapist
- Extra wide couch frame for added strength and stability
- Safe working load and lifting capacity of 320kg (705lb / 50 stone) to cater for heavy load requirements
- Push button electric operation of height and backrest angle for ease of patient positioning
- Retractable wheel system for moving and positioning of the couch
- Low patient transfer height of 42cm (approx 17")
- Breathing hole and plug for prone treatment


### Overall Dimensions

Length: 189cm (74") approx  
 Width: 102cm (40") approx  
 Height: 42cm - 91cm approx  
 (17" - 36") approx

### Upholstery Dimensions




### Options & Accessories:



Pressure Relieving Foam



Battery Back-up



Fold-down Side Support Rails



Foot Switch Operation



Paper Roll Holder

## Operational Parameters



Model 50+ Electric Version

## Standard Features

- Designed and manufactured in the UK
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Extra wide top surface of 102cm
- Electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 320kg
- Retractable castors
- Hoist clearance 140mm
- Breathing hole
- Motorised backrest

## Technical Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601 P (EN)
- Vinyl Flame Retardant to CRIB 5 (BS 7176)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]
- CE Marked
- Externally audited to ISO9001:2015 by LRQA

## Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service may be available if required
- Overseas delivery 7-10 days plus shipping time

## Upholstery Vinyl

- Flame retardant to BS 5852 Ignition Source 5 - CRIB 5 (BS 7176) Medium Hazard
- Independently tested to verify effectiveness at reducing the spread of MRSA, E.coli & other similar bacteria & fungi
- Rated at 500,000 Martindale rubs without damaging or marking the top layer, making it suitable for all high-traffic, heavy-duty areas
- Phthalate free
- Please note: If the upholstery colour is not specified, 'Almond' will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond