

CLASSIC OPERATORS CHAIRS

Used extensively in various clinical environments within the healthcare industry, these Classic Operators Chairs feature a polished aluminium base with chrome gas strut and have a tilting seat facility for additional comfort and positioning. The medium backrest has both height and angular adjustment.



● Colour Citrus Green



Model 1006 with optional fold-down armrests

● Colour Battleship Grey

Model 1005

Classic Operators Chair - Standard

Height Range 44-58cm

Seat Dimensions 48cm x 43cm (WxD)

Optional Accessories - see page 46

Model 1122 Anti-Static Castors

Model 1121 Load Locking Castors

Model 1006

Classic Operators Chairs - High Version

Height Range 57-83cm

Foot support ring fitted

Seat Dimensions 48cm x 43cm (WxD)

Model 1123 Glides for added stability

Model 1125 Manual Braked Castors