

PREMIUM MICROSCOPE CHAIRS

The new Premium Microscope Chairs are available with either a standard seat or saddle seat design. Both models feature adjustable armrests for optimum comfort and positioning when working with a microscope or similar equipment, as well as height and angle adjustable backrests and seat tilt. Model 1022, the standard seat chair, comes with armrests that can be moved forwards or backwards and are height adjustable and the armrest pads can be rotated and locked in any position. By contrast, the saddle seat chair (model 1024) features both height and width adjustable armrests with three positional settings.



● Colour Grape



● Colour Cocoa

Model 1022

Premium Microscope Chair - Standard

Height Range 45-59cm

Model 1023 Premium Microscope Chair - High

Version (not shown) - Height Range 54 - 74cm

Model 1024

Premium Microscope Saddle Chair - Standard

Height Range 47-61cm

Model 1025 Premium Microscope Saddle Chair -

High Version (not shown) - Height Range 56 - 76cm