

# PREMIUM CLINICIANS CHAIRS

A very comfortable clinical chair with injection moulded foam on the backrest and seat cushions. Both the seat and backrest sections have independent levers to set the required angle of seat tilt and the required backrest angle to maximise comfort and positioning. The foot height control option provides “hands-free” height operation which can be beneficial for certain clinical procedures.



● Colour Atlantic Blue

## **Model 1026**

---

Premium Clinicians Chair - Standard

---

Height Range 45-59cm

---

**Model 1116** Foot height control (not shown)

---

## **Model 1027**

---

Premium Clinicians Chair - High Model

---

Height Range 54-74cm

---

See page 38 for more optional accessories

See page 39 for vinyl upholstery colours