

HARRIS CLINIC CHAIRS

This high quality Clinic Chair features a contoured seat and backrest. The backrest has spring-loaded angular adjustment and is fully height adjustable in 7 different positions. The chair also provides seat tilt adjustment and there is the option of having foot height control which provides the benefit of “hands-free” height operation. A multi-positional arm support option can be added for stabilising the arms and hands during minor ops.



Model 1032 with optional foot height control

● Colour Sapphire Blue



Model 1033 with optional foot height control and multi-positional arm support

Model 1032

Premium Harris Clinic Chair - Standard

Height Range 48-64cm

Model 1033

Harris Clinic Chair - High Version

Height Range 54-68cm

Optional Accessories

Model 1117 Foot Height Control
Model 1122 Anti-Static Castors

Model 1128 Multi-Positional Arm Support
Model 1124 Foot Support Ring