

HARRIS SUPPORT CHAIRS

These chairs feature a very supportive bowl-shaped seat and sculptured backrest to support the shoulders and back. The free-float functionality of the design makes the seat and backrest follow the user's movements, therefore providing great comfort throughout lengthy seated procedures. With the single lever asynchronous mechanism the seat and back support move independently from each other and, once positioned, will both lock simultaneously when activated. The option of a foot height control provides the benefit of "hands-free" height operation.



● Colour Sapphire Blue



● Colour Beige

Model 1034

Harris Support Chair - Standard

Height Range 45-59cm

Model 1035

Harris Support Chair - High Version (not shown)

Height Range 54-74cm
