

FIXED HEIGHT PHLEBOTOMY CHAIR

This compact Fixed Height Phlebotomy Chair is ideal for use in small clinics or surgeries. The chair is supplied with a sculptured armrest to fully support the patients arm during phlebotomy procedures. The arm support can be used on either side of the chair and has both height and angular adjustment to ensure optimum positioning during these procedures.



Model 1112

Fixed Height Phlebotomy Chair with single armrest

H 800mm, Seat Height 480mm, D 420mm, W 750mm



Model 1113

Fixed Height Phlebotomy Chair with dual armrests

H 800mm, Seat Height 480mm, D 420mm, W 790mm

Options

See page 39 for vinyl upholstery colours