

SURGEONS ARM SUPPORT KIT

The Surgeons Armrest Kit can be used on most pneumatic height adjustable stools and chairs. This facility clamps onto the gas strut shaft just below the seat mechanism. The armrest poles have forward extension, height adjustment and rotation and are fitted with fully rotational forearm support pads. Once set-up, the armrests provide a stable and comfortable position for the clinician during various medical procedures.



Model 1179

Surgeons Arm Support Kit
Width between arms 50cm, Armrest extension 11cm

Forearm pads provide 360 degree rotational swivel

Armrest Height adjustment 35cm • Video available
